

FrontLine™ Mask

VIRUS ELIMINATION FACE MASK

DESIGNED FOR COMFORT
ENGINEERED TO PROTECT
REUSABLE

Custom Fabric Mask With PathoGen™ Filter

- High air flow polyester outer layer
- Cotton interior
- Easy breathing with high filtration
- Flexible interchangeable inserts
- Outer mask can be washed and sanitized

Designed for everyday use to help prevent infection spread. Advanced comfort, quick and clean filter replacement.

Engineered to capture and destroy airborne viruses, for protection beyond disposable respirators and cloth masks.

Three Disposable Filter Styles

PathoGen™ Filter

Advanced protection for higher exposure situations, vulnerable medical conditions, travel, public transportation, large crowds.

N95 Base Filter**

For general everyday use, back to school, back to office, retail and dining.

Sport Filter

For exercise, sports, outdoor activities.




High Air Flow Fabric

Easy breathe fabric accommodates replaceable filters.

Designed to fit the face, provide optimal seal and position the PathoGen™ Filter where it matters most.

Polymer stiffener directs airflow, fills out the mask creating an air chamber, holds the fabric away from the face, easy install.



LCP Medical Technologies, Inc. 7145 Boone Ave N Suite 180 Brooklyn Park, MN 55428 lcp-medical.com info@lcp-medical.com

*This product has not been evaluated by the FDA or NIOSH agencies for medical use. Patents pending. Preliminary review indicates the FrontLine™ Mask complies with ASTM F3502-21 Barrier Face Coverings Standard Level 2 when used as directed. **Use of N95 rated filter material does not imply N95 rating for mask product.

FrontLine™ Mask

lcp-medical.com


CONTACT: sales@lcp-medical.com

VIRUS ELIMINATION FACE MASK




PathoGen™ Filter

PathoGen™ Circuit embedded in N95 rated filter material



N95 Base Filter

Two layers factory tested N95 GSM filter material**




Copper Anti-Pathogen Technology

PathoGen™ Filter contains multiple surface layers of 99.9% pure elemental copper, reported in EPA Reg. No. 82012-1 to eliminate 99.9 percent of SARS-CoV-2, the virus that causes COVID-19. Airflow is directed to encounter copper surfaces during breathing, capturing pathogen-bearing moisture droplets, and disabling virus within the filter.



Sport Filter

High-flow filter material



Filter must be installed in mask for best protection. Remove filter prior to washing.



LCP Medical Technologies, Inc. 7145 Boone Ave N Suite 180 Brooklyn Park, MN 55428 lcp-medical.com info@lcp-medical.com

*This product has not been evaluated by the FDA or NIOSH agencies for medical use. Patents pending. Preliminary review indicates the FrontLine™ Mask complies with ASTM F3502-21 Barrier Face Coverings Standard Level 2 when used as directed. **Use of N95 rated filter material does not imply N95 rating for mask product.