

EasyBreathe™ Mask

CONSUMER FABRIC FACE MASK

DESIGNED FOR COMFORT
ENGINEERED TO PROTECT
REUSABLE

Custom High Filtration Fabric for Easy Breathing

- Filtration of 3 micron particles after multiple washes
- High air flow polyester outer layer
- Neoprene panels to seal air leaks
- Internal support holds fabric away from face
- Mask can be washed and sanitized

Designed for everyday use to help prevent infection spread.

Engineered with an air chamber for comfortable easy breathing without fogging glasses.


Protection beyond disposable surgical and cloth masks.



High Air Flow Fabric



Five Comfortable Sizes



Embedded nose wire, extended adjustable ear loops

No-flow neoprene lower panel enhances fit and seal for best protection



LCP Medical Technologies, Inc. 7145 Boone Ave N Suite 180 Brooklyn Park, MN 55428 lcp-medical.com info@lcp-medical.com

*This product has not been evaluated by the FDA or NIOSH agencies for medical use. Patents pending. Preliminary review indicates the EasyBreathe™ Mask complies with ASTM F3502-21 Barrier Face Coverings Standard Level 2 when used as directed. **Use of N95 rated filter material does not imply N95 rating for mask product.